OsteoBiol[®]
by TecnoSS

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- (1) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H
PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE
J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19
- (2) NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC; 10(4):264-70. EPUB 2008 JAN 30
- (3) CRESPI R, CAPPARE P, ROMANOS GE, MARIANI E, BENASCIUTTI E, GHERLONE E
CORTICOCANCELLOUS PORCINE BONE IN THE HEALING OF HUMAN EXTRACTION SOCKETS: COMBINING HISTOMORPHOMETRY WITH OSTEOBLAST GENE EXPRESSION PROFILES IN VIVO
INT J ORAL MAXILLOFAC IMPLANTS, 2011 JUL-AUG; 26(4):866-72
- (4) CASSETTA M, PERROTTI V, CALASSO S, PIATTELLI A, SINJARI B, IEZZI G
BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS
CLIN ORAL IMPLANTS RES, 2014 MAY 26 EPUB AHEAD OF PRINT
- (5) FIGUEIREDO A, COIMBRA P, CABRITA A, GUERRA F, FIGUEIREDO M
COMPARISON OF A XENOGENIC AND AN ALLOPLASTIC MATERIAL USED IN DENTAL IMPLANTS IN TERMS OF PHYSICO-CHEMICAL CHARACTERISTICS AND IN VIVO INFLAMMATORY RESPONSE
MATER SCI ENG C MATER BIOL APPL, 2013 AUG 1; 33(6):3506-13
- (6) FISCHER KR, STAVROPOULOS A, CALVO GUIRADO JL, SCHNEIDER D, FICKL S
INFLUENCE OF LOCAL ADMINISTRATION OF PAMIDRONATE ON EXTRACTION SOCKET HEALING - A HISTOMORPHOMETRIC PROOF-OF-PRINCIPLE PRE-CLINICAL IN VIVO EVALUATION
CLIN ORAL IMPLANTS RES, 2014 SEP 15 EPUB AHEAD OF PRINT
- (7) FESTA VM, ADDABBO F, LAINO L, FEMIANO F, RULLO R
PORCINE-DERIVED XENOGRAFT COMBINED WITH A SOFT CORTICAL MEMBRANE VERSUS EXTRACTION ALONE FOR IMPLANT SITE DEVELOPMENT: A CLINICAL STUDY IN HUMANS
CLIN IMPLANT DENT RELAT RES, 2011 NOV 14, EPUB AHEAD OF PRINT
- (8) CARDAROPOLI D, CARDAROPOLI G
PRESERVATION OF THE POSTEXTRACTION ALVEOLAR RIDGE: A CLINICAL AND HISTOLOGIC STUDY
INT J PERIODONTICS RESTORATIVE DENT, 2008 OCT; 28(5):469-77
- (9) BARONE A, CRESPI R, NICOLI ALDINI N, FINI M, GIARDINO R, COVANI U
MAXILLARY SINUS AUGMENTATION: HISTOLOGIC AND HISTOMORPHOMETRIC ANALYSIS
INT J ORAL MAXILLOFAC IMPLANTS, 2005, 20: 519-525
- (10) BOTTINI LP, RICCI L, PIATTELLI A, PERROTTI V, IEZZI G
BUCCO-LINGUAL CRESTAL BONE CHANGES AROUND IMPLANTS IMMEDIATELY PLACED IN FRESH SOCKETS IN ASSOCIATION OR NOT WITH PORCINE BONE: A NON-BLINDED RANDOMIZED CONTROLLED TRIAL IN HUMANS
J PERIODONTOL, 2012 OCT 29, EPUB AHEAD OF PRINT

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The Tecnos® patented manufacturing process used to obtain OsteoBiol® materials is able to achieve biocompatibility preserving part of the collagen matrix of the animal bone⁽⁶⁾ and avoiding at the same time high temperatures that would cause ceramization of the granules: the result is a unique biomaterial, consisting of mineral component and organic matrix, with a porous surface extremely similar to autogenous bone and able to resorb progressively while new bone formation takes place⁽²⁾.

mp3, a pre-hydrated cortico-cancellous bone mix with 10% collagen gel, has been developed with this innovative biotechnology and is a "ready-to-use" product.

mp3 main indication is lateral access maxillary sinus lift^(3,7,8), always in association with **Evolution** membranes, recommended to cover the antrostomy: the *mp3* syringe can be directly applied into the bony window without having to mix the *mp3* granules with saline.


Due to its collagen gel content, *mp3* allows an excellent graft stability while its hydrophilia guarantees quick blood absorption and therefore the necessary graft vascularization.

mp3 has also been successfully used in combination with **Evolution** membranes for alveolar ridge preservation⁽⁹⁾: the application of this biomaterial significantly limits the alveolar ridge width and height reduction that would naturally occur with spontaneous healing, preserving thus the alveolar ridge volume and allowing a correct second stage implant placement.

Finally, *mp3* is also indicated for horizontal augmentation (two wall defects) in combination with autogenous bone blocks⁽¹⁰⁾ or with OsteoBiol® **Lamina**⁽¹¹⁾: its cortico-cancellous composition allows a progressive resorption of

osteoclastic type, and in parallel a similar rate of new bone formation⁽²⁾.

These unique properties allow a very good graft volume preservation, a healthy new bony tissue and ultimately, a successful implant rehabilitation.

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- BARONE A, RICCI M, COVANI U, NANNMARK U, AZARMEHR I, CALVO GUIRADO JL
MAXILLARY SINUS AUGMENTATION USING PREHYDRATED CORTICO-CANCELLOUS PORCINE BONE: HISTOMORPHOMETRIC EVALUATION AFTER 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2012 JUN;14(3):373-9 EPUB 2010 MAY 11
- NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70 EPUB 2008 JAN 30
- RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATÉ SANCHEZ DE VAL JE, DELGADO RUIZ RA, NEGRI B, BARONA DORADO C
ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRAFTS USED IN MAXILLARY SINUS FLOOR ELEVATION
CLIN ORAL IMPLANTS RES, 2012 JAN 26, EPUB AHEAD OF PRINT
- CALVO GUIRADO JL, RAMIREZ FERNANDEZ MP, NEGRI B, DELGADO RUIZ RA, MATÉ SANCHEZ DE VAL JE, GOMEZ MORENO G
EXPERIMENTAL MODEL OF BONE RESPONSE TO COLLAGENIZED XENOGRAFTS OF PORCINE ORIGIN (OSTEOBIOL® mp3): A RADIOLOGICAL AND HISTOMORPHOMETRIC STUDY
CLIN IMPLANT DENT RELAT RES, 2013 FEB;15(1):143-51 EPUB 2011 MAR 31
- RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, DELGADO RUIZ RA, MATÉ SANCHEZ DE VAL JE, VICENTE ORTEGA V, MESEGUER OLMOS L
BONE RESPONSE TO HYDROXYAPATITES WITH OPEN POROSITY OF ANIMAL ORIGIN (PORCINE OSTEOBIOL® mp3) AND BOVINE (ENDOBN®): A RADIOLOGICAL AND HISTOMORPHOMETRIC STUDY
CLIN ORAL IMPLANTS RES, 2011 JUL; 22(7):767-73
- FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H
PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE
J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19
- SILVESTRI M, MARTEGANI P, D'AVENIA F, FARNETI M, CAPRI D, PAOLANTONI G, LANDI L
SIMULTANEOUS SINUS AUGMENTATION WITH IMPLANT PLACEMENT: HISTOMORPHOMETRIC COMPARISON OF TWO DIFFERENT GRAFTING MATERIALS. A MULTICENTER DOUBLE-BLIND PROSPECTIVE RANDOMIZED CONTROLLED CLINICAL TRIAL
INT J ORAL MAXILLOFAC IMPLANTS, 2013 MAR-APR;28(2):543-9
- PAGLIANI L, ANDERSSON P, LANZA M, NAPPO A, VERROCCHI D, VOLPE S, SENNERBY L
A COLLAGENATED PORCINE BONE SUBSTITUTE FOR AUGMENTATION AT NEOS IMPLANT SITES: A PROSPECTIVE 1-YEAR MULTICENTER CASE SERIES STUDY WITH HISTOLOGY
CLIN IMPLANT DENT RELAT RES, 2012 OCT;14(5):746-58 EPUB 2010 OCT 26
- BARONE A, ALDINI NN, FINI M, GIARDINO R, CALVO GUIRADO JL, COVANI U
XENOGRAFT VERSUS EXTRACTION ALONE FOR RIDGE PRESERVATION AFTER TOOTH REMOVAL: A CLINICAL AND HISTOMORPHOMETRIC STUDY
J PERIODONTOL, 2008 AUG;79(8):1370-7
- BARONE A, COVANI U
MAXILLARY ALVEOLAR RIDGE RECONSTRUCTION WITH NONVASCULARIZED AUTOGENOUS BLOCK BONE: CLINICAL RESULTS
J ORAL MAXILLOFAC SURG, 2007, 65:2039-2046
- WACHTEL H, FICKE S, HINZE M, BOLZ W, THALMAIR T
THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES
INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG;33(4):491-7

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- (1) CALVO GUIRADO JL, GOMEZ MORENO G, GUARDIA J, ORTIZ RUIZ A, PIATELLI A, BARONE A, MARTINEZ GONZALEZ JM, MESENGUER OLMO L, LOPEZ MARI L, DORADO CB
 BIOLOGICAL RESPONSE TO PORCINE XENOGRIFT IMPLANTS: AN EXPERIMENTAL STUDY IN RABBITS
 IMPLANT DENT, 2012 APR;21(2):112-7
- (2) ARCURI C, CECCHETTI F, GERMANO F, MOTTA A, SANTACROCE C
 CLINICAL AND HISTOLOGICAL STUDY OF A XENOGENIC SUBSTITUTE USED AS A FILLER IN POSTEXTRACTIVE ALVEOLUS
 MINERVA STOMATOL, 2005, 54: 351-362.
- (3) NANNMARK U, AZARMEHR I
 SHORT COMMUNICATION: COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS
 CLIN IMPLANT DENT RELAT RES, 2010 JUN 1; 12(2):161-3
- (4) BARONE A, AMERI S, COVANI U
 IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS
 EUR J IMPLANT PROSTHODONTICS, 2006, 2: 99-106
- (5) CASSETTA M, RICCI L, IEZZI G, DELL'AQUILA D, PIATELLI A, PERROTTI V
 RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN
 INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT;32(5):581-9
- (6) CALVO GUIRADO JL, PARDO ZAMORA G, SAEZ YUGUERO MR
 RIDGE SPLITTING TECHNIQUE IN ATROPHIC ANTERIOR MAXILLA WITH IMMEDIATE IMPLANTS, BONE REGENERATION AND IMMEDIATE TEMPORISATION: A CASE REPORT
 J IR DENT ASSOC, 2007 WINTER, 53(4):187-90

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- (1) BARONE A, CORNELINI R, CIAGLIA R, COVANI U. IMPLANT PLACEMENT IN FRESH EXTRACTION SOCKETS AND SIMULTANEOUS OSTEOTOME SINUS FLOOR ELEVATION: A CASE SERIES. INT J PERIODONTICS RESTORATIVE DENT, 2008 JUN; 28(3):283-9
- (2) SANTAGATA M, GUARINIELLO L, RAUSO R, TARTARO G. IMMEDIATE LOADING OF DENTAL IMPLANT AFTER SINUS FLOOR ELEVATION WITH OSTEOTOME TECHNIQUE: A CLINICAL REPORT AND PRELIMINARY RADIOGRAPHIC RESULTS. J ORAL IMPLANTOL, 2010 DEC; 36(6):485-489
- (3) COVANI U, CORNELINI R, BARONE A. BUCCAL BONE AUGMENTATION AROUND IMMEDIATE IMPLANTS WITH AND WITHOUT FLAP ELEVATION: A MODIFIED APPROACH. INT J ORAL MAXILLOFAC IMPLANTS, 2008 SEP-OCT; 23(5):841-6
- (4) CARDAROPOLI D, CARDAROPOLI G. HEALING OF GINGIVAL RECESSIONS USING A COLLAGEN MEMBRANE WITH A DEMINERALIZED XENOGRAFT: A RANDOMIZED CONTROLLED CLINICAL TRIAL. INT J PERIODONTICS RESTORATIVE DENT, 2009 FEB; 29(1):59-67
- (5) NANNMARK U, AZARMEHR I. SHORT COMMUNICATION: COLLAGENATED CORTICOCANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS. CLIN IMPLANT DENT RELAT RES, 2010 JUN 1; 12(2):161-3

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- TRUBIANI O, SCARANO A, ORSINI G, DI IORIO D, D'ARCANGELO C, PICCIRILLI M, SIGISMONDO M, CAPUTI S
 THE PERFORMANCE OF HUMAN PERIODONTAL LIGAMENT MESENCHYMAL STEM CELLS ON XENOGENIC BIOMATERIALS
 INT J IMMUNOPATHOL PHARMACOL, 2007 JAN-MAR; 20(1 SUPPL 1):87-91
- BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R
 OSTEOBIOL® INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS
 EUR J INFLAMM, 2011, VOL 9, NO. 3(S), 103-107
- ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S, PIATTELLI A
 HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL
 J PERIODONTOL 2006 DEC;77(12):1984-90
- SCARANO A, PIATTELLI A, PERROTTI V, MANZONI L, IEZZI G
 MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRICAL EVALUATION AFTER 4 AND 6 MONTHS
 CLIN IMPLANT DENT RELAT RES, 2011 MAR; 13(1):13-18
- EZZI G, DEGIDI M, PIATTELLI A, MANGANO C, SCARANO A, SHIBU JA, PERROTTI V
 COMPARATIVE HISTOLOGICAL RESULTS OF DIFFERENT BIOMATERIALS USED IN SINUS AUGMENTATION PROCEDURES: A HUMAN STUDY AT 6 MONTHS
 CLIN ORAL IMPLANTS RES, 2012 DEC;23(12):1369-76 Epub 2011 NOV 2
- BARONE A, AMERI S, COVANI U
 IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS
 EUR J IMPLANT PROSTHODONTICS, 2006, 2: 99-106
- BARONE A, TOTI P, QUARANTA A, DERCHI G, COVANI U
 THE CLINICAL OUTCOMES OF IMMEDIATE VERSUS DELAYED RESTORATION PROCEDURES ON IMMEDIATE IMPLANTS: A COMPARATIVE COHORT STUDY FOR SINGLE-TOOTH REPLACEMENT
 CLIN IMPLANT DENT RELAT RES, 2014 MAY 22 Epub ahead of print
- SCARANO A, PIATTELLI A, ASSENZA B, QUARANTA A, PERROTTI V, PIATTELLI M, IEZZI G
 PORCINE BONE USED IN SINUS AUGMENTATION PROCEDURES: A 5-YEAR RETROSPECTIVE CLINICAL EVALUATION
 J ORAL MAXILLOFAC SURG, 2010 AUG; 68(8):1869-73

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(2) FELICE P, PIANA L, CHECCHI L, CORVINO V, NANNMARK U, PIATTELLI M
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT
INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR;33(2):159-66

(3) ESPOSITO M, CANNIZZARO G, SOARDI E, PISTILLI R, PIATTELLI M, CORVINO V, FELICE P
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 6 MM-LONG, 4 MM-WIDE IMPLANTS OR BY LONGER IMPLANTS IN AUGMENTED BONE. PRELIMINARY RESULTS FROM A PILOT RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2012;5(1):19-33

(4) FELICE P, PISTILLI R, PIATTELLI M, SOARDI E, CORVINO V, ESPOSITO M
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. PRELIMINARY RESULTS FROM A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, SUMMER, 5(2):149-61

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(1) NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70

(2) CARDAROPOLI D, CARDAROPOLI G
HEALING OF GINGIVAL RECESSIONS USING A COLLAGEN MEMBRANE WITH A DEMINERALIZED XENOGRFT: A RANDOMIZED CONTROLLED CLINICAL TRIAL
INT J PERIODONTICS RESTORATIVE DENT, 2009 FEB; 29(1):59-67

(3) BARONE A, BORGIA V, COVANI U, RICCI M, PIATTELLI A, IEZZI G
FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2014 MAR 1 EPUB AHEAD OF PRINT

(4) BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U
TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING
CLIN ORAL IMPLANTS RES, 2012 JUL 12, EPUB AHEAD OF PRINT

(5) CARDAROPOLI D, CARDAROPOLI G
PRESERVATION OF THE POSTEXTRACTION ALVEOLAR RIDGE: A CLINICAL AND HISTOLOGIC STUDY
INT J PERIODONTICS RESTORATIVE DENT, 2008 OCT; 28(5):469-77

(6) RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATÉ SANCHEZ DE VAL JE, DELGADO RUIZ RA, NEGRI B, BARONA DORADO C
ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRFTS USED IN MAXILLARY SINUS FLOOR ELEVATION
CLIN ORAL IMPLANTS RES, 2012 JAN 26, EPUB AHEAD OF PRINT

(7) BARONE A, RICCI M, GRASSI RF, NANNMARK U, QUARANTA A, COVANI U
A 6-MONTH HISTOLOGICAL ANALYSIS ON MAXILLARY SINUS AUGMENTATION WITH AND WITHOUT USE OF COLLAGEN MEMBRANES OVER THE OSTEOTOMY WINDOW: RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2013 JAN;24(1):1-6. EPUB 2011 DEC 12

(8) SCARANO A, PIATTELLI A, PERROTTI V, MANZON I, IEZZI G
MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRICAL EVALUATION AFTER 4 AND 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2011 MAR; 13(1):13-18

(9) CASSETTA M, RICCI I, IEZZI G, CALASSO S, PIATTELLI A, PERROTTI V
USE OF PIEZOSURGERY DURING MAXILLARY SINUS ELEVATION: CLINICAL RESULTS OF 40 CONSECUTIVE CASES
INT J PERIODONTICS RESTORATIVE DENT, 2012 DEC;32(6):E182-8

(10) COVANI U, MARCONCINI S, CRESPI R, BARONE A
IMMEDIATE IMPLANT PLACEMENT AFTER REMOVAL OF A FAILED IMPLANT: A CLINICAL AND HISTOLOGICAL CASE REPORT
J ORAL IMPLANTOL, 2009; 35(4):189-95

(11) FELICE P, PIANA L, CHECCHI L, CORVINO V, NANNMARK U, PIATTELLI M
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT
INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR;33(2):159-66

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(3) FESTA VM, ADDABBO F, LAINO L, FEMIANO F, RULLO R
PORCINE-DERIVED XENOGRAFT COMBINED WITH A SOFT CORTICAL MEMBRANE VERSUS EXTRACTION ALONE FOR IMPLANT SITE DEVELOPMENT: A CLINICAL STUDY IN HUMANS
CLIN IMPLANT DENT RELAT RES, 2011 NOV 14, EPUB AHEAD OF PRINT

(4) WACHTEL H, FICKL S, HINZE M, BOLZ W, THALMAIR T
THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES
INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG;33(4):491-7

(5) SCARANO A, PIATTELLI M, CARINCI F, PERROTTI V
REMOVAL, AFTER 7 YEARS, OF AN IMPLANT DISPLACED INTO THE MAXILLARY SINUS. A CLINICAL AND HISTOLOGIC CASE REPORT
J OF OSSEOINTEGRATION, 2009

(6) HINZE M, VRIELINCK I, THALMAIR T, WACHTEL H, BOLZ W
ZYGOMATIC IMPLANT PLACEMENT IN CONJUNCTION WITH SINUS BONE GRAFTING: THE "EXTENDED SINUS ELEVATION TECHNIQUE". A CASE-COHORT STUDY
ORAL CRANIOFAC TISSUE ENG, 2011;1:188-197

(7) RINNA C, UNGARI C, SALTARELLI A, CASSONI A, REALE G
ORBITAL FLOOR RESTORATION
J CRANIOFAC SURG, 2005, 16(6): 968-972

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(5) FIGUEIREDO A, COIMBRA P, CABRITA A, GUERRA F, FIGUEIREDO M
COMPARISON OF A XENOGENIC AND AN ALLOPLASTIC MATERIAL USED IN DENTAL IMPLANTS IN TERMS OF PHYSICO-CHEMICAL CHARACTERISTICS AND IN VIVO INFLAMMATORY RESPONSE
MATER SCI ENG C MATER BIOL APPL, 2013 AUG 1;33(6):3506-13. EPUB 2013 MAY 3

(6) BARONE A, CRESPI R, ALDINI NN, FINI M, GIARDINO R, COVANI U
MAXILLARY SINUS AUGMENTATION: HISTOLOGIC AND HISTOMORPHOMETRIC ANALYSIS
INT J ORAL MAXILLOFAC IMPLANTS, 2005 JUL-AUG; 20(4):519-25

(7) CRESPI R, CAPPARE P, ROMANOS GE, MARIANI E, BENASCIUTTI E, GHERLONE E
CORTICOCANCELLOUS PORCINE BONE IN THE HEALING OF HUMAN EXTRACTION SOCKETS: COMBINING HISTOMORPHOMETRY WITH OSTEOBLAST GENE EXPRESSION PROFILES IN VIVO
INT J ORAL MAXILLOFAC IMPLANTS, 2011 JUL-AUG; 26(4):866-72

(8) CASSETTA M, PERROTTI V, CALASSO S, PIATTELLI A, SINJARI B, IEZZI G
BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50:50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS
CLIN ORAL IMPLANTS RES, 2014 MAY 26 EPUB AHEAD OF PRINT

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(1) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H
PHYSICOCHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE
J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19

(2) TRUBIANI O, SCARANO A, ORSINI G, DI IORIO D, D'ARCANGELO C, PICCIRILLI M, SIGISMONDO M, CAPUTI S
THE PERFORMANCE OF HUMAN PERIODONTAL LIGAMENT MESENCHYMAL STEM CELLS ON XENOGENIC BIOMATERIALS
INT J IMMUNOPATHOL PHARMACOL, 2007 JAN-MAR; 20(1 SUPPL 1):87-91

(3) NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC; 10(4):264-70. Epub 2008 JAN 30

(4) CRESPI R, CAPPARE P, ROMANOS GE, MARIANI E, BENASCIUTTI E, GHERLONE E
CORTICOCANCELLOUS PORCINE BONE IN THE HEALING OF HUMAN EXTRACTION SOCKETS: COMBINING HISTOMORPHOMETRY WITH OSTEOBLAST GENE EXPRESSION PROFILES IN VIVO
INT J ORAL MAXILLOFAC IMPLANTS, 2011 JUL-AUG; 26(4):866-72

(5) FBARONE A, RICCI M, COVANI U, NANNMARK U, AZARMEHR I, CALVO GUIRADO JL
MAXILLARY SINUS AUGMENTATION USING PREHYDRATED CORTICOCANCELLOUS PORCINE BONE: HISTOMORPHOMETRIC EVALUATION AFTER 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2012 JUN; 14(3):373-9. Epub 2010 MAY 11

(6) BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R
OSTEOBIOL® INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS
EUR J INFLAMMAT, 2011, VOL. 9, NO. 3(S), 103-107

(7) BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U
TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING
CLIN ORAL IMPLANTS RES, 2013 NOV; 24(11):1231-7. Epub 2012 JUL 12

(8) BARONE A, BORGIA V, COVANI U, RICCI M, PIATELLI A, IEZZI G
FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2014 MAR 1 Epub AHEAD OF PRINT

(9) FISCHER KR, STAVROPOULOS A, CALVO GUIRADO JL, SCHNEIDER D, FICKL S
INFLUENCE OF LOCAL ADMINISTRATION OF PAMIDRONATE ON EXTRACTION SOCKET HEALING – A HISTOMORPHOMETRIC PROOF-OF-PRINCIPLE PRE-CLINICAL IN VIVO EVALUATION
CLIN ORAL IMPLANTS RES, 2014 SEP 15 Epub AHEAD OF PRINT

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(1) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H
PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE
J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19

(2) BARONE A, RICCI M, GRASSI RF, NANNMARK U, QUARANTA A, COVANI U
A 6-MONTH HISTOLOGICAL ANALYSIS ON MAXILLARY SINUS AUGMENTATION WITH AND WITHOUT USE OF COLLAGEN MEMBRANES OVER THE OSTEOTOMY WINDOW: RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2013 JAN;24(1):1-6. EPUB 2011 DEC 12

(3) SCARANO A, PIATTELLI A, PERROTTI V, MANZON L, IEZZI G
MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRICAL EVALUATION AFTER 4 AND 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2011 MAR; 13(1):13-18

(4) CASSETTA M, PERROTTI V, CALASSO S, PIATTELLI A, SINIARI B, IEZZI G
BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS
CLIN ORAL IMPLANTS RES, 2014 MAY 26 EPUB AHEAD OF PRINT

(5) BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R
OSTEOBIOL® INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS
EUR J INFLAMMAT, 2011, VOL 9, NO. 3(S), 103-107

(6) HSU FY, CHUEH SC, WANG YJ
MICROSPHERES OF HYDROXYAPATITE/RECONSTITUTED COLLAGEN AS SUPPORTS FOR OSTEOBLAST CELL GROWTH
BIOMATERIALS 1999, 20:1931-1936

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(1) GRIFFITH LG, NAUGHTON G
TISSUE ENGINEERING-CURRENT CHALLENGES AND EXPANDING OPPORTUNITIES
SCIENCE 2002, 295:1009-14

(2) REDDI AH
MORPHOGENESIS AND TISSUE ENGINEERING OF BONE AND CARTILAGE: INDUCTIVE SIGNALS, STEM CELLS, AND BIOMIMETIC BIOMATERIALS
TISSUE ENG 2000, 6(4):351-59

(3) NAKASHIMA N, REDDI AH
THE APPLICATION OF BONE MORPHOGENETIC PROTEINS TO DENTAL TISSUE ENGINEERING
NAT BIOTECHNOL 2003, 9:1025-32

(4) SALASZNYK RM, WILLIAMS WA, BOSKEY A, BATORSKY A, PLOPPER GE
ADHESION TO VITRONECTIN AND COLLAGEN I PROMOTES OSTEOGENIC DIFFERENTIATION OF HUMAN MESENCHYMAL STEM CELLS
J BIOMED BIOTECHNOL 2004, 1:24-34

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(2) NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC; 10(4):264-70. EPUB 2008 JAN 30

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Annex III | Porcine and Equine Bone Matrix
Source: Tecnos[®] s.r.l.

Annex III | Porcine and Equine Membranes
Source: Tecnos[®] s.r.l.

Annex III | Equine Felts
Source: Tecnos[®] s.r.l.

Annex V | Porcine and Equine Bone Matrix
Source: Tecnos[®] s.r.l.

Annex V | Porcine and Equine Membranes
Source: Tecnos[®] s.r.l.

Annex V | Equine Felts
Source: Tecnos[®] s.r.l.

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2 | CASSETTA M, CALASSO S, VOZZA I, DELL'AQUILA D
REHABILITATION OF ATROPHIC ALVEOLAR CRESTS WITH CYLINDRICAL SANDBLASTED AND ACID ETCHED IMPLANTS: A PILOT STUDY
 EUR J IMPLANT PROSTHODONTICS, 2005;(3)1:133-144

3 | ARCURI C, CECCHETTI F, GERMANO F, MOTTA A, SANTACROCE C
CLINICAL AND HISTOLOGICAL STUDY OF A XENOGENIC BONE SUBSTITUTE USED AS A FILLER IN POSTEXTRACTIVE ALVEOLUS
 MINERVA STOMATOL, 2005 JUN;54(6):351-62

4 | BARONE A, CRESPI R, ALDINI NN, FINI M, GIARDINO R, COVANI U
MAXILLARY SINUS AUGMENTATION: HISTOLOGIC AND HISTOMORPHOMETRIC ANALYSIS
 INT J ORAL MAXILLOFAC IMPLANTS, 2005 JUL-AUG; 20(4):519-25

5 | RINNA C, UNGARI C, SALTARELA A, CASSONI A, REALE G
ORBITAL FLOOR RESTORATION
 J CRANIOFAC SURG, 2005 NOV; 16(6):968-72

6 | BARONE A, AMERI S, COVANI U
IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS
 EUR J IMPLANT PROSTHODONTICS, 2006,2: 99-106

7 | BARONE A, SANTINI S, SBORDONE I, CRESPI R, COVANI U
A CLINICAL STUDY OF THE OUTCOMES AND COMPLICATIONS ASSOCIATED WITH MAXILLARY SINUS AUGMENTATION
 INT J ORAL MAXILLOFAC IMPLANTS, 2006 JAN-FEB; 21(1):81-5

8 | COVANI U, BARONE A, CORNELINI R, CRESPI R
CLINICAL OUTCOME OF IMPLANTS PLACED IMMEDIATELY AFTER IMPLANT REMOVAL
 J PERIODONTOL, 2006 APR;77(4):722-7

9 | ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S, PIATTELLI A
HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL
 J PERIODONTOL, 2006 DEC;77(12):1984-90

10 | TRUBIANI O, SCARANO A, ORSINI G, DI IORIO D, D'ARCANGELO C, PICCIRILLI M, SIGISMONDO M, CAPUTI S
THE PERFORMANCE OF HUMAN PERIODONTAL LIGAMENT MESENCHYMAL STEM CELLS ON XENOGENIC BIOMATERIALS
 INT J IMMUNOPATHOL PHARMACOL, 2007 JAN-MAR; 20(1 SUPPL 1):87-91

11 | BARONE A, COVANI U
MAXILLARY ALVEOLAR RIDGE RECONSTRUCTION WITH NON- VASCULARIZED AUTOGENOUS BLOCK BONE: CLINICAL RESULTS
 J ORAL MAXILLOFAC SURG, 2007 OCT;65(10):2039-46

12 | CALVO GUIRADO JL, PARDO ZAMORA G, SAEZ YUGUERO MR
RIDGE SPLITTING TECHNIQUE IN ATROPHIC ANTERIOR MAXILLA WITH IMMEDIATE IMPLANTS, BONE REGENERATION AND IMMEDIATE TEMPORISATION: A CASE REPORT
 J IR DENT ASSOC 2007 WINTER;53(4):187-90

13 | DEL CORSO M
SOFT TISSUE RESPONSE TO PLATELET RICH FIBRIN: CLINICAL EVIDENCES
 COSMETIC DENTISTRY, 2008, 3:16-20

14 | BARONE A, SANTINI S, MARCONCINI S, GIACOMELLI L, GHERLONE E, COVANI U
OSTEOTOMY AND MEMBRANE ELEVATION DURING THE MAXILLARY SINUS AUGMENTATION PROCEDURE. A COMPARATIVE STUDY: PIEZOELECTRIC DEVICE VS. CONVENTIONAL ROTATIVE INSTRUMENTS
 CLIN ORAL IMPLANTS RES, 2008 MAY;19(5):511-5. EPUB 2008 MAR 26

15 | BARONE A, CORNELINI R, CIAGLIA R, COVANI U
IMPLANT PLACEMENT IN FRESH EXTRACTION SOCKETS AND SIMULTANEOUS OSTEOTOME SINUS FLOOR ELEVATION: A CASE SERIES
 INT J PERIODONTICS RESTORATIVE DENT, 2008 JUN; 28(3):283-9

16 | BARONE A, ALDINI NN, FINI M, GIARDINO R, CALVO GUIRADO JL, COVANI U
XENOGRAFT VERSUS EXTRACTION ALONE FOR RIDGE PRESERVATION AFTER TOOTH REMOVAL: A CLINICAL AND HISTOMORPHOMETRIC STUDY
 J PERIODONTOL, 2008 AUG;79(8):1370-7

17 | COVANI U, CORNELINI R, BARONE A
BUCCAL BONE AUGMENTATION AROUND IMMEDIATE IMPLANTS WITH AND WITHOUT FLAP ELEVATION: A MODIFIED APPROACH
 INT J ORAL MAXILLOFAC IMPLANTS, 2008 SEP-OCT; 23(5):841-6

18 | CARDAROPOLI D, CARDAROPOLI G
PRESERVATION OF THE POSTEXTRACTION ALVEOLAR RIDGE: A CLINICAL AND HISTOLOGIC STUDY
 INT J PERIODONTICS RESTORATIVE DENT, 2008 OCT; 28(5):469-77

19 | NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
 CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70. EPUB 2008 JAN 30

20 | CALVO GUIRADO JL, GOMEZ MORENO G, BARONE A, CUTANDO A, ALCARAZ BANOS M, CHIVA F, LOPEZ MARI L, GUARDIA J
MELATONIN PLUS PORCINE BONE ON DISCRETE CALCIUM DEPOSIT IMPLANT SURFACE STIMULATES OSTEOINTEGRATION IN DENTAL IMPLANTS
 J PINEAL RES, 2009, 47(2):164-72

21 | SCARANO A, PIATTELLI M, CARINCI F, PERROTTI V
REMOVAL, AFTER 7 YEARS, OF AN IMPLANT DISPLACED INTO THE MAXILLARY SINUS. A CLINICAL AND HISTOLOGIC CASE REPORT
 J OSSEOINTEGR, 2009;1(1):35-40

22 | COVANI U, MARCONCINI S, CRESPI R, BARONE A
IMMEDIATE IMPLANT PLACEMENT AFTER REMOVAL OF A FAILED IMPLANT: A CLINICAL AND HISTOLOGICAL CASE REPORT
 J ORAL IMPLANTOL, 2009; 35(4):189-95

23 | CALVO GUIRADO JL, GOMEZ MORENO G, LOPEZ MARI L, ORTIZ RUIZ AJ, GUARDIA J
ATRAUMATIC MAXILLARY SINUS ELEVATION USING THREADED BONE DILATORS FOR IMMEDIATE IMPLANTS. A THREE-YEAR CLINICAL STUDY
 MEDICINA ORAL, PATOLOGIA ORAL Y CIRUGIA BUCAL, 2010 MAR 1; 15(2):E366-70

24 | FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H
PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE
 J BIOMED MATER RES B APPL BIOMATER, 2010FEB; 92(2):409-19

25 | GRENGA PI, REALE G, COFONE C, MEDURI A, CERUTI P, GRENGA R
HESS AREA RATIO AND DIPLOPIA: EVALUATION OF 30 PATIENTS UNDERGOING SURGICAL REPAIR FOR ORBITAL BLOW-OUT FRACTURE
 OPHTHAL PLAST RECONSTR SURG, 2009 MAR-APR; 25(2):123-5

26 | CRESPI R, CAPPARE P, GHERLONE E
DENTAL IMPLANTS PLACED IN EXTRACTION SITES GRAFTED WITH DIFFERENT BONE SUBSTITUTES: RADIOGRAPHIC EVALUATION AT 24 MONTHS
 J PERIODONTOL, 2009 OCT; 80(10):1616-1621

27 | RINNA C, REALE G, FORESTA E, MUSTAZZA MC
MEDIAL ORBITAL WALL RECONSTRUCTION WITH SWINE BONE CORTEX
 J CRANIOFAC SURG, 2009 MAY; 20(3): 881-4

30 | CARDAROPOLI D, CARDAROPOLI G
HEALING OF GINGIVAL RECESSIONS USING A COLLAGEN MEMBRANE WITH A THE MINERALIZED XENOGRFT: A RANDOMIZED CONTROLLED CLINICAL TRIAL
INT J PERIODONTICS RESTORATIVE DENT, 2009 FEB; 29(1):59-67

31 | NANNMARK U, AZARMEHR I
SHORT COMMUNICATION: COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2010 JUN 1; 12(2):161-3

32 | BARONE A, RICCI M, COVANI U, NANNMARK U, AZARMEHR I, CALVO GUIRADO JL
MAXILLARY SINUS AUGMENTATION USING PRE-HYDRATED CORTICOCANCELLOUS PORCINE BONE: HYS TOMORPHOMETRIC EVALUATION AFTER 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2012 JUN; 14(3):373-9

33 | CALVO GUIRADO JL, GOMEZ MORENO G, LOPEZ MARI L, GUARDIA J, MARINEZ GONZALEZ JM, TRESGUERRES IF, PAREDES SD, FUENTES BRETO L
ACTIONS OF MELATONIN MIXED WITH COLLAGENIZED PORCINE BONE VERSUS PORCINE BONE ONLY ON OSTEOINTEGRATION OF DENTAL IMPLANTS
J PINEAL RES, 2010, 48:194-203

34 | SCARANO A, PIATTELLI A, ASSENZA B, QUARANTA A, PERROTTI V, PIATTELLI M, IEZZI G
PORCINE BONE USED IN SINUS AUGMENTATION PROCEDURES: A 5-YEAR RETROSPECTIVE CLINICAL EVALUATION
J ORAL MAXILLOFAC SURG, 2010 AUG; 68(8):1869-73

35 | ROSSI R, MORALES RS, FRASCARIA M, BENZI R, SQUADRITO N
PLANNING IMPLANTS IN THE ESTHETIC ZONE USING A NEW IMPLANT 3D NAVIGATION SYSTEM
EUR J ESTHETIC DENT, 2010 SUMMER; 5(2):172-88

37 | SCARANO A, CARINCI F, ASSENZA B, PIATTELLI M, MURMURA G, PIATTELLI A
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE USING AN INLAY TECHNIQUE WITH A XENOGRFT WITHOUT MINISCREWS AND MINIPLATES: CASE SERIES
CLIN ORAL IMPLANTS RES, 2011 OCT; 22(10):1125-30

38 | PAGLIANI L, ANDERSSON P, LANZA M, NAPPO A, VERROCCHI D, VOLPE S, SENNERBY L
A COLLAGENATED PORCINE BONE SUBSTITUTE FOR AUGMENTATION AT NEOS IMPLANT SITES: A PROSPECTIVE 1-YEAR MULTICENTER CASE SERIES STUDY WITH HISTOLOGY
CLIN IMPLANT DENT RELAT RES, 2012 OCT; 14(5):746-58

39 | SANTAGATA M, GUARINIELLO L, TARTARO G
A MODIFIED EDENTULOUS RIDGE EXPANSION (MERE) TECHNIQUE FOR IMMEDIATE PLACEMENT OF IMPLANTS. A CASE REPORT
J ORAL IMPLANTOL, 2011 MAR; 37 SPEC N°:114-9

40 | BARONE A, RICCI M, CALVO GUIRADO JL, COVANI U
BONE REMODELLING AFTER REGENERATIVE PROCEDURES AROUND IMPLANTS PLACED IN FRESH EXTRACTION SOCKETS: AN EXPERIMENTAL STUDY IN BEAGLE DOGS
CLIN ORAL IMPLANTS RES, 2011 OCT; 22(10):1131-7

41 | SCARANO A, PIATTELLI A, PERROTTI V, MANZON L, IEZZI G
MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRICAL EVALUATION AFTER 4 AND 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2011 MAR; 13(1):13-18

42 | CALVO GUIRADO JL, RAMIREZ FERNANDEZ MP, NEGRI B, DELGADO RUIZ RA, MATÉ SANCHEZ DE VAL JE, GOMEZ MORENO G
EXPERIMENTAL MODEL OF BONE RESPONSE TO COLLAGENIZED XENOGRFTS OF PORCINE ORIGIN (OSTEOBIOL® MP3): A RADIOLOGICAL AND HISTOMORPHOMETRIC STUDY
CLIN IMPLANT DENT RELAT RES, 2013 FEB; 15(1):143-51

43 | RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, DELGADO RUIZ RA, MATÉ SANCHEZ DE VAL JE, VICENTE ORTEGA V, MESEGUER OLMO S
BONE RESPONSE TO HYDROXYAPATITES WITH OPEN POROSITY OF ANIMAL ORIGIN (PORCINE OSTEOBIOL® MP3) AND BOVINE (ENDOBON®): A RADIOLOGICAL AND HISTOMORPHOMETRIC STUDY
CLIN ORAL IMPLANTS RES, 2011 JUL; 22(7):767-73

44 | CRESPI R, CAPPARÈ P, ROMANOS GE, MARIANI E, BENASCIOTTI E, GHERLIONE E
CORTICOCANCELLOUS PORCINE BONE IN THE HEALING OF HUMAN EXTRACTION SOCKETS: COMBINING HISTOMORPHOMETRY WITH OSTEOBLAST GENE EXPRESSION PROFILES IN VIVO
INT J ORAL MAXILLOFAC IMPLANTS, 2011 JUL-AUG; 26(4):866-72

45 | HINZE M, VRIELINCK L, THALMAIR T, WACHTEL H, BOLZ W
ZYGOMATIC IMPLANT PLACEMENT IN CONJUNCTION WITH SINUS BONE GRAFTING: THE "EXTENDED SINUS ELEVATION TECHNIQUE". A CASE-COHORT STUDY
ORAL CRANIOFAC TISSUE ENG, 2011; 1:188-197

46 | IEZZI G, DEGIDI M, PIATTELLI A, MANGANO C, SCARANO A, SHIBLI JA, PERROTTI V
COMPARATIVE HISTOLOGICAL RESULTS OF DIFFERENT BIOMATERIALS USED IN SINUS AUGMENTATION PROCEDURES: A HUMAN STUDY AT 6 MONTHS
CLIN ORAL IMPLANTS RES, 2012 DEC; 23(12):1369-76


47 | SLOTTE C, LINDFORS N, NANNMARK U
SURGICAL RECONSTRUCTION OF PERI-IMPLANT BONE DEFECTS WITH PREHYDRATED AND COLLAGENATED PORCINE BONE AND COLLAGEN BARRIERS: CASE PRESENTATIONS
CLIN IMPLANT DENT RELAT RES, 2013 OCT; 15(5):714-23

48 | BARONE A, RICCI M, GRASSI RF, NANNMARK U, QUARANTA A, COVANI U
A 6-MONTH HISTOLOGICAL ANALYSIS ON MAXILLARY SINUS AUGMENTATION WITH AND WITHOUT USE OF COLLAGEN MEMBRANES OVER THE OSTEOTOMY WINDOW: RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2013 JAN; 24(1):1-6

49 | SANTAGATA M, GUARINIELLO L, RAUSO R, TARTARO G
IMMEDIATE LOADING OF DENTAL IMPLANT AFTER SINUS FLOOR ELEVATION WITH OSTEOTOME TECHNIQUE: A CLINICAL REPORT AND PRELIMINARY RADIOGRAPHIC RESULTS
J ORAL IMPLANTOL, 2010 DEC; 36(6):485-489

50 | FESTA VM, ADDABBO F, IAINO L, FEMIANO F, RULLO R
PORCINE-DERIVED XENOGRFT COMBINED WITH A SOFT CORTICAL MEMBRANE VERSUS EXTRACTION ALONE FOR IMPLANT SITE DEVELOPMENT: A CLINICAL STUDY IN HUMANS
CLIN IMPLANT DENT RELAT RES, 2013 OCT; 15(5):707-13

51 | RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATÉ SANCHEZ DE VAL JE, DELGADO RUIZ RA, NEGRI B, BARONA DORADO C
ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRFTS USED IN MAXILLARY SINUS FLOOR ELEVATION
CLIN ORAL IMPLANTS RES, 2013 MAY; 24(5):523-30

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52 | CASSETTA M, RICCI L, IEZZI G, DELL'AQUILA D, PIATTELLI A, PERROTTI V
RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN
INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT; 32(5):581-9

53 | BARONE A, ORLANDO B, CINGANO L, MARCONCINI S, DERCHI G, COVANI U
A RANDOMIZED CLINICAL TRIAL TO EVALUATE AND COMPARE IMPLANTS PLACED IN AUGMENTED VS. NON-AUGMENTED EXTRACTION SOCKETS. A 3-YEAR EVALUATION
J PERIODONTOL, 2012 JUL; 83(7):836-46

54 | CALVO GUIRADO JL, MATÉ SANCHEZ JE, DELGADO RUIZ L, RAMIREZ FERNANDEZ MP
CALCULATION OF BONE GRAFT VOLUME USING 3D RECONSTRUCTION SYSTEM
MED ORAL PATOL CIR BUCAL, 2011 MAR 1, 16(2):E260-4

55 | ESPOSITO M, CANNIZZARO G, SOARDI E, PISTILLI R, PIATTELLI M, CORVINO V, FELICE P
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 6 MM-LONG, 4 MM-WIDE IMPLANTS OR BY LONGER IMPLANTS IN AUGMENTED BONE. PRELIMINARY RESULTS FROM A PILOT RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2012 SPRING; 5(1):19-33

56 | FELICE P, PIANA L, CHECCHI L, PISTILLI R, PELLEGRINO G
VERTICAL RIDGE AUGMENTATION OF THE ATROPHIC POSTERIOR MANDIBLE WITH A 2-STAGE INLAY TECHNIQUE: A CASE REPORT
IMPLANT DENT, 2012 JUN; 21(3):190-5

57 | BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U
TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING
CLIN ORAL IMPLANTS RES, 2013 NOV; 24(11):1231-7

58 | CASSETTA M, RICCI L, IEZZI G, CALASSO S, PIATTELLI A, PERROTTI V
USE OF PIEZOSURGERY DURING MAXILLARY SINUS ELEVATION: CLINICAL RESULTS OF 40 CONSECUTIVE CASES
INT J PERIODONTICS RESTORATIVE DENT, 2012 DEC; 32(6):E182-8

59 | BOTTINI LP, RICCI L, PIATTELLI A, PERROTTI V, IEZZI G
BUCCO-LINGUAL CRESTAL BONE CHANGES AROUND IMPLANTS IMMEDIATELY PLACED IN FRESH SOCKETS IN ASSOCIATION OR NOT WITH PORCINE BONE: A NON-BLINDED RANDOMIZED CONTROLLED TRIAL IN HUMANS
J PERIODONTOL, 2012 OCT 29, EPUB AHEAD OF PRINT

60 | BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R
OSTEOBIOL® INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS
EUR J INFLAMMAT, 2011, VOL. 9, NO. 3 (S), 103-107

61 | FELICE P, PIANA L, CHECCHI L, CORVINO V, NANNMARK U, PIATTELLI M
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT
INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR; 33(2):159-66

62 | FICKL S, JOCKEL-SCHNEIDER Y, LINCKE T, BECHTOLD M, FISCHER KR, SCHLAGENHAUF U
PORCINE DERMAL MATRIX FOR COVERING OF RECESSION TYPE DEFECTS: A CASE SERIES
QUINTESENCE INT, 2013; 44(3):243-6

63 | SILVESTRI M, MARTEGANI P, D'AVENIA F, FARNETI M, CAPRI D, PAOLANTONI G, LANDI L
SIMULTANEOUS SINUS AUGMENTATION WITH IMPLANT PLACEMENT: HISTOMORPHOMETRIC COMPARISON OF TWO DIFFERENT GRAFTING MATERIALS. A MULTICENTER DOUBLE-BLIND PROSPECTIVE RANDOMIZED CONTROLLED CLINICAL TRIAL
INT J ORAL MAXILLOFAC IMPLANTS, 2013 MAR-APR; 28(2):543-9

64 | WACHTEL H, FICKL S, HINZE M, BOLZ W, THALMAIR T
THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES
INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG; 33(4):491-7

65 | RODRIGUEZ JG, ELDBANY RM
VERTICAL SPLITTING OF THE MANDIBULAR BODY AS AN ALTERNATIVE TO INFERIOR ALVEOLAR NERVE LATERALIZATION
INT J ORAL MAXILLOFAC SURG, 2013 SEP; 42(9):1060-6

66 | FIGUEIREDO A, COIMBRA P, CABRITA A, GUERRA F, FIGUEIREDO M
COMPARISON OF A XENOGENEIC AND AN ALLOPLASTIC

MATERIAL USED IN DENTAL IMPLANTS IN TERMS OF PHYSICO-CHEMICAL CHARACTERISTICS AND IN VIVO INFLAMMATORY RESPONSE
MATER SCI ENG C, MATER BIOL APP, 2013 AUG 1; 33(6):3506-13

67 | FELICE P, PISTILLI R, PIATTELLI M, SOARDI E, CORVINO V, ESPOSITO M
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. PRELIMINARY RESULTS FROM A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, SUMMER; 5(2):149-61

68 | TRAINI T, PIATTELLI A, CAPUTI S, DEGIDI M, MANGANO C, SCARANO A, PERROTTI V, IEZZI G
REGENERATION OF HUMAN BONE USING DIFFERENT BONE SUBSTITUTE BIOMATERIALS
CLIN IMPLANT DENT RELAT RES, 2013 MAY 17, EPUB AHEAD OF PRINT

69 | CALVO GUIRADO JL, GOMEZ MORENO G, GUARDIA J, ORTIZ RUIZ A, PIATTELLI A, BARONE A, MARTÍNEZ GONZÁLEZ JM, MESENGUER OLMO L, LÓPEZ MARI L, DORADO CB
BIOLOGICAL RESPONSE TO PORCINE XENOGRFT IMPLANTS: AN EXPERIMENTAL STUDY IN RABBITS
IMPLANT DENT, 2012 APR; 21(2):112-7

70 | KOLMAS J, SZWAJA M, KOŁODZIEJSKI W
SOLID-STATE NMR AND IR CHARACTERIZATION OF COMMERCIAL XENOGENEIC BIOMATERIALS USED AS BONE SUBSTITUTES
J PHARM BIOMED ANAL, 2012 MAR 5; 61:136-41

71 | PISTILLI R, FELICE P, PIATTELLI M, GESSAROLI M, SOARDI E, BARALUSSE C, BUTI J, CORVINO V, ESPOSITO M
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. ONE-YEAR RESULTS FROM A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2013 WINTER; 6(4):343-357

72 | FICKL S, NANNMARK U, SCHLAGENHAUF U, HÜRZELER M, KEBSCHULL M
PORCINE DERMAL MATRIX IN THE TREATMENT OF DEHISCENCE-TYPE DEFECTS – AN EXPERIMENTAL SPLIT-MOUTH ANIMAL TRIAL
CLIN ORAL IMPLANTS RES, 2014 FEB 19 EPUB AHEAD OF PRINT

73 | BARONE A, BORGIA V, COVANI U, RICCI M, PIATTELLI A, IEZZI G
FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2014 MAR 1 EPUB AHEAD OF PRINT

74 | BARONE A, RICCI M, ROMANOS GE, TONELLI P, ALFONSI F, COVANI U
BUCCAL BONE DEFICIENCY IN FRESH EXTRACTION SOCKETS: A PROSPECTIVE SINGLE COHORT STUDY
CLIN ORAL IMPLANTS RES, 2014 MAR 31 EPUB AHEAD OF PRINT

75 | BARONE A, TOTI P, QUARANTA A, DERCHI G, COVANI U
THE CLINICAL OUTCOMES OF IMMEDIATE VERSUS DELAYED RESTORATION PROCEDURES ON IMMEDIATE IMPLANTS: A COMPARATIVE COHORT STUDY FOR SINGLE-TOOTH REPLACEMENT
CLIN IMPLANT DENT RELAT RES, 2014 MAY 22 EPUB AHEAD OF PRINT

76 | CASSETTA M, PERROTTI V, CALASSO S, PIATTELLI A, SINJARI B, IEZZI G
BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS
CLIN ORAL IMPLANTS RES, 2014 MAY 26 EPUB AHEAD OF PRINT

77 | GHENO E, PALERMO A, BUFFOLI B, RODELLA LF
THE EFFECTIVENESS OF THE USE OF XENOGENEIC BONE BLOCKS MIXED WITH AUTOLOGOUS CONCENTRATED GROWTH FACTORS (CGF) IN BONE REGENERATION TECHNIQUES
J OSSEointegration 2014; 6(2):37-42

78 | FAUSI G, GALLI M, VITTORINI-VELASQUEZ P, GALLEGOS-RIVERA JC, MINASI R, DE BIASE A, DI PAOLO C
USE OF 3D CARTILAGE SCAFFOLDS FOR THE STABILIZATION OF IMPLANTS AND BONE REGENERATION WITH THE FIT-LOCK TECHNIQUE
ACTA ODONTOLATINOAM 2013; 26(3):167-172

79 | FISCHER KR, FICKL S, MARDAS N, BOZEC L, DONOS N
STAGE-TWO SURGERY USING COLLAGEN SOFT TISSUE GRAFTS: CLINICAL CASES AND ULTRASTRUCTURAL ANALYSIS
QUINTESENCE INT, 2014 NOV-DEC; 45(10):853-60


80 | FISCHER KR, STAVROPOULOS A, CALVO GUIRADO JL, SCHNEIDER D, FICKL S
INFLUENCE OF LOCAL ADMINISTRATION OF PAMIDRONATE ON EXTRACTION SOCKET HEALING – A HISTOMORPHOMETRIC PROOF-OF-PRINCIPLE PRE-CLINICAL IN VIVO EVALUATION
CLIN ORAL IMPLANTS RES, 2014 SEP 15 EPUB AHEAD OF PRINT

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
| МАТЕРИАЛ | УПАКОВКА | ОПИСАНИЕ | РАЗМЕРЫ И ВЕС | СВИНОЕ ПРОИСХОЖДЕНИЕ | КОНСКОЕ ПРОИСХОЖДЕНИЕ |
|-------------------------------|------------|--------------------|----------------------|----------------------|-----------------------|
| BONE SUBSTITUTES | | | | | |
| Gen-Os | 1 Vial | DRIED GRANULES | 0.25 g | M1052FS | M1052FE |
| Gen-Os | 1 Vial | DRIED GRANULES | 0.5 g | M1005FS | M1005FE |
| Gen-Os | 1 Vial | DRIED GRANULES | 1.0 g | M1010FS | M1010FE |
| Gen-Os | 1 Vial | DRIED GRANULES | 2.0 g | M1020FS | M1020FE |
| mp3 | 1 Syringe | BONE MIX | 1.0 cc | A3005FS | A3005FE |
| mp3 | 3 Syringes | BONE MIX | 3x0.25 cc (0.75 cc) | A3075FS | |
| mp3 | 3 Syringes | BONE MIX | 3x0.5 cc (1.5 cc) | A3015FS | A3015FE |
| mp3 | 3 Syringes | BONE MIX | 3x1.0 cc (3.0 cc) | A3030FS | A3030FE |
| Putty | 1 Syringe | BONE PASTE | 0.5 cc | HPT09S | HPT09E |
| Putty | 1 Syringe | BONE PASTE | 1.0 cc | HPT61S | HPT61E |
| Putty | 3 Syringes | BONE PASTE | 3x0.25 cc (0.75 cc) | HPT32S | HPT32E |
| Putty | 3 Syringes | BONE PASTE | 3x0.5 cc (1.5 cc) | HPT35S | HPT35E |
| Gel 40 | 1 Syringe | BONE GEL | 0.5 cc | 05GEL40S | 05GEL40E |
| Gel 40 | 3 Syringes | BONE GEL | 3x0.5 cc (1.5 cc) | 15GEL40S | 15GEL40E |
| MEMBRANES AND BARRIERS | | | | | |
| Evolution | 1 Blister | DRIED / FINE | 20x20x (0.3) mm | | EVO2LLE |
| Evolution | 1 Blister | DRIED / FINE | 30x30x (0.3) mm | | EVO3LLE |
| Evolution | 1 Blister | DRIED / FINE | Oval 25x35x (0.3) mm | | EVOLLE |
| Evolution | 1 Blister | DRIED / STANDARD | 20x20x (0.5) mm | EM02HS | EVO2HHE |
| Evolution | 1 Blister | DRIED / STANDARD | 30x30x (0.5) mm | EM03HS | EVO3HHE |
| Evolution | 1 Blister | DRIED / STANDARD | Oval 25x35x (0.5) mm | EM00HS | |
| Soft Cortical Lamina | 1 Blister | DRIED / FINE | 25x25x (0.5) mm | LS25FS | LS25FE |
| Soft Cortical Lamina | 1 Blister | DRIED / FINE | 35x25x (0.5) mm | LS23FS | LS23FE |
| Soft Cortical Lamina | 1 Blister | DRIED / STANDARD | 30x30x (3.0) mm | LS03SS | |
| Curved Lamina | 1 Blister | DRIED / MEDIUM | 35x35x (1.0) mm | LS10HS | LS10HE |
| SPECIFIC PRODUCTS | | | | | |
| Apatos Mix | 1 Vial | DRIED GRANULES | 0.5 g | A1005FS | A1005FE |
| Apatos Mix | 1 Vial | DRIED GRANULES | 1.0 g | A1010FS | A1010FE |
| Apatos Mix | 1 Vial | DRIED GRANULES | 2.0 g | A1020FS | A1020FE |
| Apatos Cortical | 1 Vial | DRIED GRANULES | 0.5 g | AC1005FS | |
| Apatos Cortical | 1 Vial | DRIED GRANULES | 1.0 g | AC1010FS | |
| Sp-Block | 1 Blister | DRIED BLOCK / NORM | 10x10x10 mm | | BN0E |
| Sp-Block | 1 Blister | DRIED BLOCK / NORM | 10x10x20 mm | | BN1E |
| Sp-Block | 1 Blister | DRIED BLOCK / NORM | 10x20x20 mm | | BN2E |
| Sp-Block | 1 Blister | DRIED BLOCK / NORM | 35x10x5 mm | | BN8E |
| Dual-Block CURVED | 1 Blister | DRIED BLOCK / SOFT | 20x15x5 mm | STS7S | |
| Dual-Block CURVED | 1 Blister | DRIED BLOCK / NORM | 20x10x5 mm | STN5S | |
| Special | 1 Blister | DRIED / X-FINE | 20x20x (0.2) mm | EM02LS | EM02LE |
| Special | 1 Blister | DRIED / X-FINE | 30x30x (0.2) mm | EM03LS | EM03LE |
| Duo-Teck | 1 Blister | DRIED | 20x20x (0.2) mm | | DT020 |
| Duo-Teck | 6 Blister | DRIED FELT | 25x25x (0.15) mm | | DTN625 |
| Derma | 1 Blister | DRIED / STANDARD | 30x30x (2.0) mm | ED03SS | |
| Derma | 1 Blister | DRIED / FINE | 25x25x (0.9) mm | ED25FS | |
| Derma | 1 Blister | DRIED / STANDARD | 7x5x (2.0) mm | ED75SS | |
| Derma | 1 Blister | DRIED / STANDARD | 15x5x (2.0) mm | ED15SS | |

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